

# Dress / Fashion Orthotics — TL - Blue/2100™ Orthotics

Graphite orthotics are best suited for patients requiring maximum control with less bulk than traditional sport devices. We recommend the TL-2100 material for your patients involved in high stress/low impact activities. These products are available in several flexibilities based on patient weight, age and activity.

## Fashion Step I

- Polypropylene shell (1/8" only)
- Intrinsic forefoot posting
- Intrinsic rearfoot posting
- Fits dress flats / pumps up to 2.5" (exact heel height required)



## Fashion Step II

- TL-2100 / Blue carbon graphite shell
- Intrinsic forefoot posting
- Intrinsic rearfoot posting
- Fits dress flats / pumps up to 2.5" (exact heel height required)

### Common pathologies:

- Mild forefoot varus / valgus
- Mild rearfoot valgus
- Plantar fasciitis
- Tarsal coalition



## Execu-Step

- Polypropylene shell (5/32" only)
- Intrinsic forefoot posting
- Extrinsic rearfoot posting
- Fits laced oxfords and most loafers

### Common pathologies:

- Mild to moderate forefoot varus / valgus
- Plantar fasciitis
- Heel spurs
- Tarsal coalition

## Eliminator 2100

- TL - Blue/2100 graphite shell
- Intrinsic forefoot posting
- Intrinsic rearfoot posting
- Fits athletic shoes, full cut laced oxfords and most flats / loafers

### Common pathologies:

- Mild forefoot varus / valgus
- Mild rearfoot valgus
- Plantar fasciitis
- Heel spurs
- Tarsal coalition

## Eliminator Plus

- TL - Blue/2100 carbon graphite shell
- Intrinsic forefoot posting
- Extrinsic rearfoot posting - 65 durometer SBR rubber
- Fits athletic shoes, full cut laced oxfords, some flats and most loafers.

### Common pathologies:

- Mild to moderate forefoot varus / valgus
- Moderate to severe rearfoot varus / valgus
- Plantar fasciitis
- Heel spurs
- Tarsal coalition

JSB does not recommend dress orthotics for heels over 2" (There is an additional charge for dress orthotics over 2")

### ELIMINATOR 2100 and PLUS Material breakdown / patient weight

- 1.8mm Semi-Flexible TL-2100.....Up to 120#
- 2.2mm Semi-Rigid TL-2100.....121# to 185#
- 2.8mm Rigid TL-2100.....186# to 250#
- 3.2mm Ultra TL2100.....251# and up

### FASHION STEP II Material breakdown / patient weight

- 1.8mm Semi-Flexible TL-2100.....Up to 150#
- 2.2mm Semi-Rigid TL-2100.....151# to 200#
- 2.8mm Rigid TL-2100.....201# to 250#
- Ultra TL-2100.....251# and up