

ELIMINATOR LINE - SPORT ORTHOTICS

(POLYPROPYLENE/HIGH DENSITY POLYETHYLENE)

The Eliminator line is our group of all-purpose sport orthotics. JSB's system of varying material thickness based on patient's weight, activity and physiology, assures your patient both maximum orthotic performance and optimum fit. Please be sure to complete the area on the prescription form asking for weight, height, age, sex, occupation and shoe data.



Eliminator I

- Polypropylene shell
- Intrinsic forefoot posting
- Intrinsic rearfoot posting
- Fits athletic shoes and most flat dress shoes

Common pathologies:

- Mild forefoot varus / valgus
- Mild rearfoot valgus / valgus
- Morton's Neuroma

Eliminator II

- Polypropylene shell
- Intrinsic forefoot posting
- Extrinsic rearfoot posting - 65 durometer SBR rubber
- Fits athletic shoes and full cut laced oxfords

Common pathologies:

- Mild to moderate forefoot varus / valgus
- Moderate to severe rearfoot varus / valgus
- Calcaneal apophysitis

- Plantar fasciitis
- Heel spurs
- Tarsal coalition

Eliminator III

The Eliminator III incorporates an extrinsic forefoot post for maximum biomechanical control

- Polypropylene shell
- Extrinsic forefoot posting - 65 durometer SBR rubber
- Extrinsic rearfoot posting - 65 durometer SBR rubber
- Fits athletic shoes and full cut laced oxfords

Common pathologies:

- Moderate to severe forefoot varus / valgus
- Moderate to severe rearfoot varus / valgus
- Lisfranc fracture / syndrome

Eliminator IV

- High density polyethylene shell
- Intrinsic forefoot posting
- Extrinsic rearfoot posting 65 durometer SBR rubber
- Soft Evazote arch fill / over 250# Firm Pelite arch fills
- 1mm EVA bottom covers
- Fits athletic shoes and full cut laced oxfords

Common pathologies:

- Mild to moderate forefoot varus / valgus
- Moderate to severe rearfoot varus / valgus
- Plantar fasciitis
- Heel Spurs
- Tarsal Coalition

Eliminator I and II
Material breakdown / patient weight

- 3.2mm Polypro (1/8").....Up to 150#
- 4.0mm Polypro (5/32").....151# to 225#
- 4.7mm Polypro (3/16").....226# to 300#
- 6.0mm Polypro (1/4").....301# and Up

Eliminator III
Material breakdown / patient weight

- 3.2mm Polypro (1/8").....Up to 175#
- 4.0mm Polypro (5/32").....176# to 225#
- 4.7mm Polypro (3/16").....226 to 300#
- 6.0mm Polypro (1/4").....301# and Up

Eliminator IV
Material breakdown / patient weight

- 3.2 HDPE (1/8").....All Weights